

CLAP facial roller, 3.5



Lyapko roller applicator, 40x51 mm

Rating: Not Rated Yet

Price

29,00 €

Discount:

[Ask a question about this product](#)

Manufacturer: [ARGO company](#)

Description Elastic roller, with needles made of different metals that are necessary for the body: zinc, copper, iron, nickel, silver.
Step of the needles: 3.5 mm, size: 40x51 mm

It is recommended to use in cosmetology, it is most convenient to influence face and head zones for pain and spasm relief, (for trigeminal neuralgia, facial neuritis treatment).

Due to its size, it allows you to handle the space around the damaged areas (trauma, post-surgical scars, fractures), varicose veins (without exacerbation).

"The needle rain" with facial roller recovers tired nervous system, relaxes facial muscles, rejuvenates your skin. Use the facial roller, facing to a mirror, - and you will by all means want to smile to yourselves.

Application: The first 5 minutes of relatively uncomfortable stinging sensation from the effects of the applicator turn into a comfortable feeling with the emergence of a powerful pleasant heat, "vibration", tingling. Subsequently, it may be a sense of total relaxation and sleepiness, sometimes passing into a healthy, restful sleep.

- applicator exposure time for analgesic, tonic effect, to enhance performance (daytime) and debilitated patients and children - 5-15 minutes;
- to express the analgesic and hypnotic effect (evening) - for 20-25 minutes. Course duration is 1-3 weeks, and sometimes - 2-3 days. Repeated course is possible after 1-3 weeks.

WARNING: after discomfort for 10-15 minutes, applicator should be removed and used only after 5-10 hours or the next day, capturing a smaller area. Applicator should cause only a pleasant experience! Needle of roller in contact with the skin don't damage it. Applicator can be used at any age.

Important: please [read this article for details](#), watch the videos and learn indications and contraindications for use of applicators!

Reviews

There are yet no reviews for this product.